# Investigating the role of insect vector glycosylation in African sleeping sickness transmission: Characterisation of procainamide-labelled tsetse fly saliva N-glycans Ludger Radoslaw P. Kozak<sup>1</sup>, Katherine Wongtrakul-Kish<sup>1</sup>, Christopher Williams<sup>2</sup>, Samirah Perally<sup>2</sup>, Claire Rose<sup>2</sup>, Daniel I.R. Spencer<sup>1</sup>, and Alvaro Acosta-Serrano<sup>2</sup>

<sup>1</sup> Ludger Ltd, Culham Science Centre, Oxfordshire, UK; <sup>2</sup> Liverpool School of Tropical Medicine, Liverpool L3 5QA, UK

# Introduction

African trypanosomiasis or sleeping sickness occurs in sub-Saharan African countries and is transmitted through the saliva of the haematophagus insect vector (tsetse fly) during feeding. The causative agents of this disease are trypanosome parasites of the species *Trypanosoma brucei*. Although sustained efforts to curb infection have resulted in a decrease of cases from 9,878 in 2009 to 7,216 in 2012<sup>1</sup>, efforts to identify new drug targets to treat or prevent infection continue. Salivary glycoproteins have been reported to facilitate host infection through binding and transport of vector-borne diseases to host tissues, and also participate in host responses such as inflammation and immune response<sup>2</sup>. This role in infection presents

2. Comparison of procainamide labelled N-glycans from naïve and T. brucei-infected tsetse fly saliva glycoproteins by HILIC-SPE ESI-MS/MS



new opportunities to identify the key mediators of transmission as well as to increase our understanding of

the role of salivary glycans in haematophagus insects. The analysis of fly salivary glycoproteins however is challenging due to small sample volumes despite collection from several hundred flies, together with the need for high sensitivity.

# Aims

1. To determine a suitable glycomics workflow coupling HPLC with mass spectrometry to combat the problem of small sample amounts in tsetse fly saliva. Specifically, to compare and evaluate the use of two different labels (2-AB and procainamide) for released N-glycans from a standard glycoprotein (IgG) to ensure highly sensitive and reliable mass spectrometric and fluorescence detection

2. To investigate the role of N-glycosylation in vector based trypanosome infection by characterising the *N*-glycome of *T. brucei*-infected and naïve tsetse salivary glycoproteins

# Methods

**TSETSE FLY SALIVA:** Naïve vs. T. brucei-infected samples



#### **Bovine Fetuin standard glycoprotein:** 2-AB vs. Procainamide labelling of released *N*-glycans

•Pauci mannose and high mannose-type structures showed high expression the N-glycomes of both naïve and

*T. brucei*-infected tsetse fly saliva glycoproteins with (Man)3(GlcNAc)2-Proc the dominant structure in both

•Low abundance hybrid-type glycans were also detected in both samples

(not visible in average MS profile)



Tsetse fly<sup>3</sup> Trypanosoma brucei<sup>4</sup> •Naïve (~5µg/µl) and *T. brucei*-infected (~1µg/µl) tsetse fly saliva collection •In-gel release of *N*-glycans with PNGase F [E-PNG01 from Ludger]

Labelled with either: Released N-glycan labelling via reductive amination with Procainamide 2-Aminobenzamide (2-AB) Labelled with: using LudgerTag<sup>™</sup> 2-AB Glycan Glycan Labelling Kit *[LT-*Labelling Kit [LT-KAB-VP24] *KPROC-VP24*] •Released N-glycan labelling via reductive amination with Procainamide using Ludger Procainamide Glycan Labelling Kit •2-AB labelled *N*-glycan clean up using LudgerClean<sup>™</sup> T1 glycan clean up using • Procainamide labelled *N*-glycan clean up using LudgerClean<sup>™</sup> S Cartridges [LC-T1-A6] LudgerClean<sup>™</sup> S Cartridges [LC-S-A6] Analysed Analysed Analysed using: using: using: •HILIC-UHPLC-FLR coupled online with positive ion mode ESI-MS/MS and HILIC-SPE ESI-MS using a Bruker amaZon speed ion trap BRUKER •Manual data analysis GF

•50 µg fetuin was reduced and alkylated with DTT and IAA respectively

• *N*-glycans were released in-solution enzymatically with PNGase F [E-PNG01] using a Hamilton Robotics StarLet high-throughput liquid-handling platform

> • Released *N*-glycan labelling via reductive amination with using Ludger Procainamide • Procainamide labelled N-

MS/MS fragmentation contains diagnostic ions that aid in N-glycan structural characterisation

#### 3. Total structures identified in tsetse saliva glycoproteins by MS/MS fragmentation

HPLC Peak Id	Glucose <sup>-</sup> Units .	Tsetse Fly Saliva N-glycans Identified by HILIC-UHLPC ESI-MS/MS								HPLC		Т	Tsetse Fly Saliva N-glycans Identified by HILIC-UHLPC ESI-MS/MS (cont.)						
		Composition			Theoretical $[m/z]^{1+}$	Theoretical [ <i>m/z</i> ] <sup>2+</sup>	Detected [ <i>m/z</i> ] <sup>1+</sup>	Detected $[m/z]^{2+}$	Proposed Structure	Peak Id	Glucose Units	Composition		Theoretical	Theoretical	Detected	Detected	Proposed	
		Hex HexNAc Fuc		Hex								HexNAc	: Fuc	$[m/z]^{1+}$	$[m/z]^{2+}$	$[m/z]^{1+}$	$[m/z]^{2+}$	Structure	
1	3.21	2	2	0	968.46	484.73	968.47	n.d.	PROC	9	6 87	6	2	0	1616 67	808 84	1616 61	<u> </u>	-FROC
2	4.17	3	2	0	1130.51	565.76	1130.49	565.74			0.07		Z						
3	4.62	3	2	1	1276.57	638.79	1276.52	638.77	нос	10	7.79	7	2	0	1778.72	889.86	1778.68	889.84	-PROC
4	4.76	3	3	0	1333.59	667.30	1333.57	667.29											
5	5.00	4	2	0	1292.56	646.78	1292.54	646.74	-	11	8.53	8	2	0	1940.77	970.89	n.d.	970.87	
6	5.54	4	3	0	1495.64	748.32	1495.65	748.30		12	8.66	8	2	0	1940.77	970.89	n.d.	970.87	x2
7	6.00	5	2	0	1454.61	727.81	1454.57	727.79		13	9.35	9	2	0	2102.83	1051.92	n.d.	1051.90	
8	6.46	5	3	0	1657.69	829.35	1657.64	829.33		n.d. =	not dete	cted							

# Conclusions

• Here we present a procainamide labelling system suitable for HPLC-MS based glycomic analysis of insect saliva *N*-glycans small sample amounts showing high sensitivity

### Results

[LT-KPROC-VP24]

Cartridges [LC-S-A6]

#### **1.** Comparison of 2-AB and procainamide labelled bovine fetuin released *N*-glycans: **HILIC-UHPLC-ESI-MS:** HILIC-UHPLC-FLR:



Procainamide and 2-AB showed good comparability however procainamide labelling displayed 2x higher

#### fluorescence and 40x higher MS signal intensity compared to 2-AB labelled N-glycans.

• These results suggest that upon colonisation and maturation of trypanosomes in the tsetse salivary

glands, there are no detectable changes in the glycosylation of tsetse salivary glycoproteins.

• The presence of high levels of mannosylated structures may influence the half-life in blood of

### pharmacologically active salivary components.

#### Acknowledgements

This work was partially supported by funding from the Wellcome Trust and GlycoPar EU FP7 Marie Curie Initial Training Network

#### References

- 1. http://www.who.int/mediacentre/factsheets/fs259/en/
- 2. R. G. Titus, J. V. Bishop, and J. S. Mejia, "The immunomodulatory factors of arthropod saliva and the potential for these factors to serve as vaccine targets to prevent pathogen transmission," Parasite Immunology, vol. 28, no. 4, pp. 131–141, 2006
- 3. http://www.raywilsonbirdphotography.co.uk/Galleries/Invertebrates/vectors/Tsetse\_Fly.html
- 4. http://www.york.ac.uk/cii/media-library/